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



c) Navigator > Acquire at Items, Initial Actions after Moving Stage 选择 Run script“map”, Primary Task 选择 “Acquire map image or montage”, 点击“Go”开始拍中地图。

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Tilt Angle Specifications

Tilt to End at Base increment

Increment proportional to cosine of angle Total # of tilts: 35

Run series in two directions from Return with Walk Up

Magnification for anchor: um Use View for anchor

Use dose-symmetric scheme

Delay time after tilting by basic increment: seconds

Change some settings as a function of angle

Size, Image Shift, and Recentering Specifications

Magnification: Binning: Pixel size: 0.207 nm

At start and when needed, track at magnification:

Limit image shift to um at start

Recenter beam from Low Dose tracking images

Auto-center beam every minutes

Auto-center beam in second half of series at degrees

Refine the zero-loss peak every minutes

Beam Intensity or Exposure Control

Change Record exposure time instead of beam intensity

Keep exposure time constant

Vary exposure to be 2 x higher at 60 deg (1/cos to 1/1 power)

Set exposure to keep mean counts of Record images at

The image in buffer A implies a Record zero-tilt mean of 1518

Taper counts down to above +/- degrees

Keep exposure below current value (use if exposure already set up)

Do not increase exposure above value for first saved image

Use current exposure at zero tilt (increase to start series)

Total dose: 35.0 x dose at 0 deg; 35.0 x dose at -51 deg

Change other exposures proportional to Record change

Change drift settling proportional to exposure change

Camera and Acquisition Parameters

Ceta Falcon 4EC-EF

Average Record dark references

If taking frames, get initial partial sum of frames

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Do autofocus at least every degrees Skip Autofocusing

Focus every time above degrees

Repeat autofocus if focus differs from prediction by > um

Do autofocus when error in focus prediction exceeds um

Check autofocus, stop if measured < X actual defocus

Stop if autofocus exceeds absolute focus limits

Stop if autofocus change exceeds um

Initial and Final Actions

Align to image now in Buffer A

Refine eucentricity

Leave mid-tilt anchor near plus or minus degrees

Use mid-tilt anchor that you have left in buffer T

Close column valves at end of series

Tracking Control Parameters

Wait for drift to settle

Repeat Record if percentage of field lost is more than

Stop and ask before repeating Record in Low Dose mode

Stop if Autoalign shifts more than % of image size

Get tracking image when error in X/Y prediction is > % of field

Track: Before After Before & after Autofocusing

Stop after aligning tracking image to allow manual adjustment

Align only with tracking images, not Record images

Align with Preview before getting new track reference

Get new track reference if Record alignment differs by > %

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